WorkFlow.md 6/22/2020

## Still in development

## **Imaging Workflow**

- 1. Turn on camera
- 2. Wait for EOS to open
- 3. Check EOS settings
  - 1. Confirm that renameimage. exe is 'registered' to EOS
    - 1. Preferences > Linked Software > Register > Browse select renameimage.exe
    - 2. Full path -

C:Users\Image\Documents\PNWHerbaria\Scripts\renameimage\renameimag
e.exe

- 4. Select portal folder
  - 1. In EOS Preferences > Destination Folder
  - 2. Select one of the portal folders C: Users\Image\Desktop\Imaging\
  - 3. Personal projects and other special collections that are not going on Symbiota should be imaged in the **Random** folder if you want the images backed up to the server (Sassafrass). Any folders in **Random** will be moved as is to the local backup drive, not sorted into portal folders.
- 5. Take photo
  - 1. A small box should pop up with the current file name and an place to enter the new name
- 6. Scan barcode
  - 1. Barcode should appear as new file name
  - 2. Click the "rename" button to rename the file
  - 3. The small window will close after the file is renamed.

## Repeat 5. & 6. for all photos

- 7. Batch edit photos in DPP
- 8. Convert from CR2 to JPG
  - 1. Now you are done, make sure to record on paper what specimens you imaged
- 9. Upload CSV
  - 1. CSV files ready for upload are on in the CSVLogs/ folder. On workstation 2 at LSU.
  - 2. After uploading CSV files please move them to the CSVLogged/ folder.
  - 3. If file names are edited and re-uploaded a new CSV file will be made for the date the original file was uploaded. This CSV can be combined with or replace the older CSV in the CSVLogged folder.

WorkFlow.md 6/22/2020

## Behind the scenes

Task Scheduler set to run following scripts.

RunOrganize.sh - runs organizeIncomingImages.py at 8PM every day.

organizeIncomingImages.py - checks barcodes in incoming folders, moves to appropriate folder(based on barcode) on storage drive. Put files with names that are too long or too short in BadBarcode folder.

WakeUp.sh - wakes up the computer to sit and wait for 5 minutes at 10PM, this allows the server to connect to the storage drive, in order to copy over files.

rsyncDaily.sh - runs on server to copy files to server and create derivatives. Then it calls LogtoCSV.py to turn logs into csv files to upload to portals.

LogtoCSV.py - takes all rsync logs that were modified in the last 24 hours and creates or adds to csv files that have not been uploaded to a portal yet.